

TS/ 252/11

ROSS WR LIQUID

Created on: Feb 2013

Washroom Cleaner and Sanitiser Concentrate

CONSTITUTION Aqueous mixture of non-corrosive ingredients.

APPEARANCE Crimson Red liquid

SPECIFIC GRAVITY AT 25°C 0.997 ± 0.100

pH- VALUE 6.00-8.00

MISCIBILITY Miscible in cold water

- HIGHLIGHTS**
- Economical in use – product being concentrated makes 100 ltrs user solution per litre of **ROSS WR LIQ**
 - Product is abrasive -free and corrosive-free, does not cause scratch marks on surfaces that are cleaned.
 - Efficiently removes oil stains, dirt, and water marks.
 - Ideal for cleaning bathroom fittings.
 - One step cleaning and sanitisation.
 - Safe for use on marble and granite.
 - Quickly removes oils, stains, dirt and water-marks.

GUIDELINES

At 10-20/ liter concentration, **Ross WR LIQUID** will clean and sanitise normal surface.

At 20-40/ liter concentration, **Ross WR LIQUID** will clean and sanitise heavy soiled surface.

Apply **Ross WR LIQUID** to all surfaces to be cleaned. Scrub gently to remove dirt, stains etc.

The flush handle, the toilet seat cover and the other bowl need to be cleaned and wiped with **Ross WR LIQUID** after the toilet bowl has been cleaned.

For bathroom floor clean and mop using **Ross WR LIQUID**.

For daily cleaning **Ross WR LIQUID** can be used with spray and

The information and data contained herein has been compiled based on information we believe reliable. Users should thoroughly test all applications and independently conclude satisfactory performance before commercialization, as these recommendations are non-binding.

Users assume all liabilities for use of the chemicals. We are not liable for any advice which we may have failed to give.
 *TM- TRADEMARKS OF ROSSARI BIOTECH INDIA P. LTD., ROSSARI ENZYMES, A DIVISION OF ROSSARI BIOTECH INDIA P. LTD.;
 NEUTRON IMPEX P. LTD., A SUBSIDIARY OF ROSSARI BIOTECH INDIA P. LTD.*

ROSSARI BIOTECH LIMITED

201 A & B, Ackruti Corporate Park, LBS Marg, Next to GE Gardens, Kanjurmarg (West), Mumbai - 400 078, Tel: +91-22-61233800 / 25777604 / 605, Fax: +91-22-61233846 / 871 / 25796982
 E-mail: info@rossarimail.com, Website: www.rossari.com

wipe method.

STORAGE

Ross WR LIQUID has to be stored in an air-tight, cool & dry place, away from direct heat and sunlight

SAFE HANDLING

Do not mix with any other chemicals other than advised by company

PACKING

Ross WR LIQUID is available in 2 x 5 ltrs packs.

The information and data contained herein has been compiled based on information we believe reliable. Users should thoroughly test all applications and independently conclude satisfactory performance before commercialization, as these recommendations are non-binding.

Users assume all liabilities for use of the chemicals. We are not liable for any advice which we may have failed to give.

"TM- TRADEMARKS OF ROSSARI BIOTECH INDIA P. LTD., ROSSARI ENZYMES, A DIVISION OF ROSSARI BIOTECH INDIA P. LTD.; NEUTRON IMPEX P. LTD., A SUBSIDIARY OF ROSSARI BIOTECH INDIA P. LTD."

ROSSARI BIOTECH LIMITED

201 A & B, Ackruti Corporate Park, LBS Marg, Next to GE Gardens, Kanjurmarg (West), Mumbai - 400 078, Tel: +91-22-61233800 / 25777604 / 605, Fax: +91-22-61233846 / 871 / 25796982
E-mail: info@rossarimail.com, Website: www.rossari.com